

SIGVARIS
GROUP



2024 Product Catalog

Graduated Compression
Socks & Hosiery for
Veno-Lymphatic Disorders

sigvaris

5 reasons patients love Sigvaris products

- 1. Great Fit:** Our PFS2 sizing system was designed to fit the broadest range of patients, providing accurate compression and exceptional comfort.
- 2. Easy to Put On:** Sigvaris Medical hosiery glides on easier than other compression brands because we use double-covered spandex inlay yarns.
- 3. Easy to Launder:** Products are machine washable and machine dryable. No need to air dry!
- 4. Six Month Guarantee:** SIGVARIS GROUP guarantees the compression of our Medical products for six months, even when machine washed and dried.
- 5. Color and Material Choices:** More choices of colors, patterns and even premium natural fibers like Sea Island Cotton, Linen, and Merino Wool than any other medical compression brand.





Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our passion for helping people feel their best is why we've been a global market leader for 160 years and counting.

Our products restore, assist, and enhance people's lives.

We are renowned for making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

For more information, please visit:

sigvaris.com



Live your life. Be yourself.

	Product	Page	Styles	mmHg
Essential	Cotton 230	6	<ul style="list-style-type: none"> ● Calf, Calf w/grip-top, Thigh w/grip-top ○ Calf, Thigh w/grip-top 	20–30, 30–40
	Opaque 860 <i>for women</i>	8	<ul style="list-style-type: none"> ● Calf, Petite Calf, Calf Plus, Calf w/grip-top, Thigh w/grip-top, Petite Thigh w/grip-top, Extra-Large Thigh w/grip-top, Pantyhose, Petite Pantyhose, Plus-Sized Pantyhose ○ Calf, Calf Plus, Calf w/grip-top, Thigh w/grip-top, Extra-Large Thigh w/grip-top, Thigh w/waist attachment 	20–30, 30–40
	Opaque 860 <i>for men</i>	10	<ul style="list-style-type: none"> ● Calf, Calf Plus, Calf w/grip-top, Thigh w/grip-top, Extra-Large Thigh w/grip-top ○ Calf, Calf Plus, Calf w/grip-top, Thigh w/grip-top, Extra-Large Thigh w/grip-top, Thigh w/waist attachment 	20–30, 30–40
Style	Sea Island Cotton 220	12	<ul style="list-style-type: none"> ● Calf 	20–30
	Linen 250	14	<ul style="list-style-type: none"> ● Calf 	20–30
	Medium Sheer 750 <i>for women</i>	16	<ul style="list-style-type: none"> ● Calf, Thigh w/grip-top, Pantyhose 	20–30
	Sheer 780 <i>for women</i>	18	<ul style="list-style-type: none"> ● Calf, Thigh w/grip-top, Pantyhose 	15–20, 20–30, 30–40
	Microfiber 820 <i>for men</i>	20	<ul style="list-style-type: none"> ● Calf, Calf w/grip-top, Thigh w/grip-top, Extended Foot Sizes 	15–20, 20–30, 30–40
	Microfiber Patterns 830	22	<ul style="list-style-type: none"> ● Calf 	20–30
Motion	Soft Opaque 840 <i>for women</i>	24	<ul style="list-style-type: none"> ● Calf, Calf w/grip-top, Thigh w/grip-top, Pantyhose, Maternity Pantyhose ○ Calf, Thigh w/grip-top, Pantyhose 	15–20, 20–30, 30–40
	Cushioned Cotton 360	26	<ul style="list-style-type: none"> ● Calf 	20–30
	High Tech 412	28	<ul style="list-style-type: none"> ● Calf 	20–30
	Thermoregulating Wool 422	30	<ul style="list-style-type: none"> ● Calf 	20–30
	Natural Rubber 500	32	<ul style="list-style-type: none"> ○ Calf, Thigh, Thigh w/grip-top, Thigh w/waist attachment, Pantyhose 	30–40, 40–50, 50–60
	Secure 550	34	<ul style="list-style-type: none"> ● Calf w/grip-top, Thigh w/grip-top ○ Calf w/grip-top, Thigh w/grip-top 	20–30, 30–40, 40–50
Specialty	Secure Armsleeve, Glove, and Gauntlet 560	36	<ul style="list-style-type: none"> ○ Armsleeve w/grip-top (Glove and Gauntlet available separately) 	15–20, 20–30
	Secure Lite Armsleeve 570	38	<ul style="list-style-type: none"> ○ Armsleeve w/grip-top (Glove and Gauntlet available separately) 	15–20, 20–30
	Advance Armsleeve 910	40	<ul style="list-style-type: none"> ● Armsleeve w/grip-top w/Gauntlet 	15–20, 20–30
	Diabetic Compression Socks 602	42	<ul style="list-style-type: none"> ● Calf 	18–25

To order:

Call: 1.800.322.7744, ext. 1
Fax: 1.800.481.5488

Online Dealer Portal:
sigvaris-online.com

Product		Page	Styles	mmHg
Medically Complex Edema	Compreflex® Standard Calf & Foot	46	● Lower extremity inelastic wrap	20–50
	Compreflex Standard Calf	46	● Lower extremity inelastic wrap	20–50
	Compreflex Transition Calf	48	● Lower extremity inelastic wrap	20–50
	Compreflex Complete Calf	48	● Lower extremity inelastic wrap	20–50
	Compreflex Contain Calf	50	● Lower extremity inelastic wrap	20–50
	Coolflex® Standard Calf	52	● Lower extremity inelastic wrap	20–50
	Assorted Liners	53	● Lower extremity	
	Assorted Foot Garments	54	● Foot	
Dynaven	Opaque for Women & Men	56	● Calf, Calf w/grip-top, Thigh w/grip-top, Pantyhose ○ Calf, Calf w/grip-top, Thigh w/grip-top, Pantyhose	15–20, 20–30, 30–40
	Opaque for Men	56	● Calf, Calf w/grip-top	15–20, 20–30, 30–40
	Sheer for Women	58	● Calf, Thigh w/grip-top, Pantyhose ○ Calf, Thigh w/grip-top	15–20, 20–30
	Cushioned for Women & Men	58	● Calf	15–20, 20–30
Well Being	Casual Cotton 146 / 186	60	● Calf	15–20
	Microfiber Shades 143 / 183	61	● Calf	15–20
	Sea Island Cotton 151/ 191	62	● Calf	15–20
	All-Season Merino Wool 152 / 192	63	● Calf	15–20
	Sheer Fashion 120	64	● Calf, Thigh w/grip-top, Pantyhose ○ Calf	15–20
	Soft Silhouette Leggings 170	66	○ Legging	15–20
	Athletic Recovery Socks 401	67	● Calf	15–20
	Performance Sleeves 412V	68	○ Calf Sleeve	20–30
	Merino Outdoor Socks 421	69	● Calf	15–20
	Eversoft Diabetic Sock 160	70	● Calf	8–15

Accessories
see page 72

References
see page 76

For socks and hosiery:

● = Closed Toe

○ = Open Toe

Measuring Instructions
see page 78

Ordering Instructions
see page 79

Essential

COTTON 230

Essential Cotton combines a classic ribbed design and all-day cotton comfort that is ideal for sensitive skin.

Compression levels

20–30mmHg • 30–40mmHg

Features and benefits

- Supima® cotton for highest level of comfort and breathability
- Ideal for sensitive skin and allergies
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fiber content


- 50% Nylon, 25% Spandex, 25% Supima Cotton (content may vary slightly by style)
- Non-latex

Best suited for

- Containment of various types of edema including venous edema, phlebolymphe'dema, and early stages of lymphedema
- Ideal for dress or casual shoes
- Ideal for long-distance travel



232NMLM66
20–30mmHg, Closed Toe Thigh
Medium Long in Light Beige



Silicone dotted grip-top band
available for calf and thigh

Essential Cotton 230 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"-8.5" 18-21.5cm		8.5"-10" 21.5-25.5cm		10"-11.5" 25.5-29cm		11.5"-14" 29-36cm	
Calf circumference	11"-15.5" 28-39cm		12.5"-17.5" 31.5-45cm		14"-20" 35.5-51cm		16"-24" 41-61cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"-26" 45.5-66cm		19"-28.5" 48-73cm		20.5"-33" 52-84cm		NA	
Leg length: Women & Open Toe	Short	Long	Short	Long	Short	Long		
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm		
Leg length: Men	Short	Long	Short	Long	Short	Long		
	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm		
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12 Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

Styles

Colors





Closed Toe = ●
Open Toe = ○

White
00

Navy
10

Light Beige
66

Black
99

				
Calf (C)	●	●	●○	●
Calf w/grip-top (C/S)** Silicone band grip-top			●*	●**
Thigh w/grip-top (N) Sensinnov band grip-top			●○	●

* 20-30mmHg only

** 30-40mmHg calf with grip-top only available in Black, 30-40mmHg extra large calf only available in Black with grip-top

Product code example: 232CMLW00

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
23	2	C	ML	W	00

Product code example: 232CLSW66/S

Series	mmHg Level	Style	Size	Gender / Open Toe	Color	Grip-top
23	2	C	LS	W	66	/S

Essential

OPAQUE 860

Essential Opaque offers our broadest range of sizes and styles, making it your go-to choice for everyday hosiery.

(For Men's Essential Opaque, see page 10)

Compression levels

20–30mmHg • 30–40mmHg

Features and benefits

- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles
- The Pantyhose Plus style is a great option during pregnancy

Fiber content


- 64% Nylon, 36% Spandex (content may vary slightly by style)
- Non-latex

Best suited for

- Containment of various types of edema including venous edema, phlebolympheidema, and early stages of lymphedema
- Post-vein treatment including thermal and non-thermal techniques
- Patients looking for petite-sized and plus-sized options



862CSLW99/S
20–30mmHg, Closed Toe Calf w/grip-top
Small Long in Black



Silicone dotted grip-top band
available for calf and thigh

Women's Essential Opaque 860 sizing chart *(for Men's Sizing Chart, see page 11)*

	Small			Medium			Large		X Large	
Ankle circumference	7"-8.5" 18-21.5cm			8.5"-10" 21.5-25.5cm			10"-11.5" 25.5-29cm		11.5"-14" 29-36cm	
Calf circumference	11"-15.5" 28-39cm			12.5"-17.5" 31.5-45cm			14"-20" 35.5-51cm		16"-24" 41-61cm	
Calf Circumference: Calf Plus w/grip-top	NA			15.5"-21" 39.5-53cm			18"-24" 45.5-61cm		20"-28" 51-71cm	
Calf length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤14" 35.5cm	≤16" 41cm	>16" 41cm	≤14" 35.5cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"-26" 45.5-66cm			19"-28.5" 48-73cm			20.5"-33" 52-84cm		25"-37.5" 64-95cm	
Leg length: Women & Open Toe	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤25.5" 65cm	≤28.5" 72cm	>28.5" 72cm	≤25.5" 65cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Hip circumference guidelines	32.5"-44.5" 83-112.5cm			34.5"-50" 87-127cm			38.5"-59.5" 98-150.5cm		NA	
Size	SP	SS	SL	MP	MS	ML	LS	LL	XS	XL
Calf Plus w/grip-top					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Women's Shoe Sizes: Petite Leg 4-8, Short Leg 5.5-9.5, Long Leg 8-12

Plus (Maternity) Hip Circumference Guidelines: Small ≤56" (142cm), Medium ≤62" (157.5cm), Large ≤70" (178cm)

Women's Styles *(for Men's, see page 11)*

Colors

Closed Toe = ●
Open Toe = ○

Honey
33

Golden
36

Light
Beige
66

Black
99

	Honey 33	Golden 36	Light Beige 66	Black 99
Women's Calf (C)	●	●	○	●●
Women's Calf Petite Sizes (C) Only available in 20-30mmHg, not available with silicone band grip-top	●	●		●
Women's Calf w/grip-top (C/S) Silicone band grip-top			●●	●●
Women's Calf Plus w/grip-top (C/P) Only available with silicone band grip-top			●●	●●
Women's Thigh w/grip-top (N) Silicone band grip-top*	●	●	○	●●
Women's Thigh w/grip-top Petite (N)** Only available in 20-30mmHg	●	●		●
Women's Pantyhose (P) Not available in XS or XL sizes	●	●		●
Women's Pantyhose Petite (P)** Only available in 20-30mmHg	●	●		●
Women's Pantyhose Plus (P/P) Not available in XS or XL sizes			●	●
Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)			○	

* Women's Thigh w/grip-top (N) and Open Toe Thigh w/grip-top (N) sizes XS and XL are only available in Light Beige and Black

** Thigh w/Waist Attachment not available in XS or XL sizes. Petite sizes not available in L and XL sizes.

Product code example: 862NSLW36

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
86	2	N	SL	W	36

Product code example: 862PLSW66/P

Series	mmHg Level	Style	Size	Gender / Open Toe	Color	Plus size
86	2	P	LS	W	66	/P

Essential

OPAQUE 860

Essential Opaque offers durability you can trust in our most extensive size range, making it your go-to choice for everyday socks and stockings.

(For Women's Essential Opaque, see page 8)

Compression levels

20–30mmHg • 30–40mmHg

Features and benefits

- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Men's calf and thigh products feature ribbed fabric; open toe products are a smooth hosiery fabric
- Flat, low friction toe seam

Fiber content

- 64% Nylon, 36% Spandex (content may vary slightly by style)
- Non-latex

Best suited for

- Containment of various types of edema including venous edema, phlebolympheidema, and early stages of lymphedema
- Post-vein treatment including thermal and non-thermal techniques
- Patients looking for plus-sized options



862CMLM99
20–30mmHg, Closed Toe Calf
Medium Long in Black



Silicone dotted grip-top band
available for calf and thigh

Men's Essential Opaque 860 sizing chart *(for Women's Sizing Chart, see page 9)*

	Small		Medium		Large		X Large	
Ankle circumference	7"-8.5" 18-21.5cm		8.5"-10" 21.5-25.5cm		10"-11.5" 25.5-29cm		11.5"-14" 29-36cm	
Calf circumference	11"-15.5" 28-39cm		12.5"-17.5" 31.5-45cm		14"-20" 35.5-51cm		16"-24" 41-61cm	
Calf circumference: Calf Plus w/grip-top	NA		15.5"-21" 39.5-53cm		18"-24" 45.5-61cm		20"-28" 51-71cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"-26" 45.5-66cm		19"-28.5" 48-73cm		20.5"-33" 52-84cm		25"-37.5" 64-95cm	
Leg length: Men	Short	Long	Short	Long	Short	Long	Short	Long
	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm
Leg length: Open Toe	Short	Long	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Size	SS	SL	MS	ML	LS	LL	XS	XL
Calf Plus w/grip-top			MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Men's shoe sizes: Short leg 7-11, Long leg 9.5-14



Men's Styles *(for Women's, see page 9)*

Colors

Closed Toe = ●
Open Toe = ○

Light
Beige
66

Black
99

		
Men's Calf (C)	●○	●○
Men's Calf w/grip-top (C/S) Silicone band grip-top	●○	●○
Men's Calf Plus w/grip-top (C/P) Only available with silicone band grip-top	●○	●○
Men's Thigh w/grip-top (N) Silicone band grip-top	●○	●○
Thigh w/Waist Attachment (W)* Indicate leg LEFT (/L) or RIGHT (/R)	○	

*Thigh w/Waist Attachment not available in XS or XL sizes.

Product code example: 862NLL099

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
86	2	N	LL	O	99

Product code example: 863CSSM66/S

Series	mmHg Level	Style	Size	Gender / Open Toe	Color	Grip-top
86	3	C	SS	M	66	/S

Style

SEA ISLAND COTTON 220

Style Sea Island Cotton socks are made with the world's rarest cotton, ensuring a luxurious feel that is as soft as cashmere.

Compression levels

20–30mmHg (also available in 15–20mmHg, see p. 62)

Features and benefits

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fiber content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Patients who prefer natural fibers for everyday wear
- Long-distance travel



222CSLW10
20–30mmHg, Closed Toe Calf
Small Long in Navy



222CLLW99
20–30mmHg, Closed Toe Calf
Large Long in Black

Style Sea Island Cotton 220 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm		11.5"–14" 29–36cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm		16"–24" 41–61cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Style	Colors	
Closed Toe = ●	Navy 10	Black 99
		
Calf (C)	●	●

Product code example: 222CMSM10					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
22	2	C	MS	M	10

Style

LINEN 250

Our stylish Linen socks in four gorgeous colors will make your legs look and feel great!

Compression levels

20–30mmHg

Features and benefits

- Lightweight and breathable for all day comfort
- Natural fibers are ideal for sensitive skin
- Crafted from fibers of the flax plant, our linen line is sustainable and stylish

Fiber content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Patients who prefer natural fibers for everyday wear
- Patients who prefer fashionable color options
- Long-distance travel



252CMLM205
20–30mmHg, Closed Toe Calf
Medium Long in Light Grey



252CSLW201
20–30mmHg, Closed Toe Calf
Small Long in Denim

Style Linen 250 sizing chart

	Small		Medium		Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm	
Calf length	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Style	Colors			
Closed Toe = ●	Brown 111	Denim 201	Light Grey 205	Lavender 207
				
Women's Calf (C)		●	●	●
Men's Calf (C)	●		●	

Product code example: 252CSLW207

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
25	2	C	SL	W	207

Style

MEDIUM SHEER 750

Thanks to the durability of Style Medium Sheer, you will have the confidence to wear sheer hosiery every day.

Compression levels

20–30mmHg

Features and benefits

- Fashionable, durable, medium sheer hosiery for everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles

Fiber content

- 68% Nylon, 32% Spandex (content may vary slightly by style)
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment
- Daily wear for business or casual
- Active professionals needing a durable sheer



752CLSW33
20–30mmHg, Closed Toe Calf
Large Short in Honey



Silicone dotted grip-top band
available for thigh

Style Medium Sheer 750 sizing chart

	Small		Medium		Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm	
Calf length	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"–26" 45.5–66cm		19"–28.5" 48–73cm		20.5"–33" 52–84cm	
Leg length	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Hip circumference guidelines	32.5"–44.5" 83–112.5cm		34.5"–50" 87–127cm		38.5"–59.5" 98–150.5cm	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Styles	Colors		
Closed Toe = ●	<div><div>Honey33</div><div>Nightshade94</div><div>Black99</div></div>		
	<div><div></div><div></div><div></div></div>		
Calf (C)	●	●	●
Thigh w/grip-top (N) Silicone band grip-top	●	●	●
Pantyhose (P)	●	●	●

Product code example: 752NLLW99					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
75	2	N	LL	W	99

Product code example: 752PSLW33					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
75	2	P	SL	W	33

Style

SHEER 780

Style Sheer hosiery offers fashionable sheerhness you can feel comfortable in, making it suitable for every occasion.

Compression levels

15–20mmHg • 20–30mmHg
30–40mmHg

Features and benefits

- Fashionable for professional and dress wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles
- Available in a diverse range of skin tones

Fiber content

- 20–30mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30–40mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment
- Patients who prefer fashionable options
- Ideal for business or fashion



781CLSW34
15–20mmHg, Closed Toe Calf
Large Short in Warm Sand



Sensinno band grip-top
available for thigh

Style Sheer 780 sizing chart

	Small		Medium		Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm	
Calf length	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"–26" 45.5–66cm		19"–28.5" 48–73cm		20.5"–33" 52–84cm	
Leg length	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Hip circumference guidelines	32.5"–44.5" 83–112.5cm		34.5"–50" 87–127cm		38.5"–59.5" 98–150.5cm	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Styles	Colors							
Closed Toe = ● Open Toe = ○	Honey 33	Warm Sand 34*	Golden 36	Toasted Almond 37*	Café 73	Mocha 85	Dark Navy 08	Black 99
								
Calf (C)	●○	●○	●○	●○	●○	●○	●	●○
Thigh w/grip-top (N) Sensinnov band grip-top	●○	●○	●○	●○	●○	●○	●	●○
Pantyhose (P)	●○	●○	●○	●○	●○	●○	●	●○

* 15-20mmHg & 20-30mmHg ONLY

Product code example: 782PMLW37					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
78	2	P	ML	W	37

Product code example: 783NSSW85					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
78	3	N	SS	W	85

Style

MICROFIBER 820

Style Microfiber socks offer microfiber softness for your professional and everyday wear.

Compression levels

15–20mmHg • 20–30mmHg
30–40mmHg

Features and benefits


- 100% microfiber against the skin for an ultra-soft feel
- Excellent breathability and moisture management for added comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fiber content

- 87% Nylon, 13% Spandex (Calf)
- 68% Nylon, 32% Spandex (Thigh)
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Patients who prefer fashionable options
- Extended foot sizes for men
- Long-distance travel



Silicone dotted grip-top band
available for calf



Sensinnov grip-top band
available for thigh

Style Microfiber 820 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm		11.5"–14" 29–36cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm		16"–24" 41–61cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"–26" 45.5–66cm		19"–28.5" 48–73cm		20.5"–33" 52–84cm		N/A	
Leg length	Short	Long	Short	Long	Short	Long		
	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm	≤30" 76cm	>30" 76cm		
Size	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
Extended Foot Size				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: Short leg 7–11, Long leg 9.5–14 (for extended foot sizes: Long leg 12.5–17.5)

* MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30mmHg only.

Extended foot sizes not available with grip top

Styles




Colors

Closed Toe = ●

Tan Khaki
32

Steel Grey
93***

Black
99

			
Calf (C)	●	●	●
Calf w/grip-top (C/S)** Silicone band grip-top	●		●
Thigh w/grip-top (N)** Sensinnov band grip-top	●		●

** Calf w/grip-top and Thigh w/grip-top not available in 15–20mmHg

*** Steel Grey color only available in 15–20mmHg and 20–30mmHg

Product code example: 822CMLM32

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
82	2	C	ML	M	32

Style

MICROFIBER PATTERNS 830

Style Microfiber Patterns socks combine ultimate comfort with a variety of fashionable colors and designs to spice up your wardrobe.

Compression levels

20–30mmHg (also available in 15–20mmHg, see p. 61)

Features and benefits

- 100% microfiber against the skin for an ultra-soft feel
- Fashionable, durable socks for everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off

Fiber content

- 87% Nylon, 13% Spandex (content can vary slightly by style)
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Patients who prefer fashionable options
- Long-distance travel



832CMSW02
20–30mmHg, Closed Toe Calf
Medium Short in Mariner Stripe
















832CMLM91
20–30mmHg, Closed Toe Calf
Medium Long in Graphite Chevron

Style Microfiber Patterns 830 sizing chart

	Small		Medium		Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm	
Calf length	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Style	Colors					
Closed Toe = ●	Mariner Stripe 02	Cranberry Chevron 05	Pink Stripe 15	Graphite Stripe 17	Onyx Stripe 19	Graphite Heather 27
						
Women's Calf (C)	●	●	●		●	●
Men's Calf (C)				●	●	●

Style	Colors						
Closed Toe = ●	Black Plaid 42	Purple Argyle 46	Graphite Argyle 47	Royal Blue Argyle 48	Wisteria Windows 61	Espresso Windows 89	Graphite Chevron 91
							
Women's Calf (C)		●	●	●	●		●
Men's Calf (C)	●		●	●		●	●

Product code example: 832CSSW05					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
83	2	C	SS	W	05

Product code example: 832CMSM91					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
83	2	C	MS	M	91

Style

SOFT OPAQUE 840

Style Soft Opaque hosiery blends a fashionable opaque look with long-lasting microfiber softness for your professional and everyday wear.

Compression levels

15–20mmHg • 20–30mmHg
30–40mmHg

Features and benefits

- Fashionable opaque hosiery for professional and everyday wear
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- Cotton-lined gusset for pantyhose styles

Fiber content

- 74% Nylon, 26% Spandex (content can vary slightly by style)
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- Pre- and post-vein procedure
- Patients who prefer fashionable options
- Pregnancy
- Long distance travel



842CMLW60
20–30mmHg, Closed Toe Calf
Medium Long in Mulberry










Sensinnov grip-top band
available for thigh

Style Soft Opaque 840 sizing chart

	Small		Medium		Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm	
Calf length	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"–26" 45.5–66cm		19"–28.5" 48–73cm		20.5"–33" 52–84cm	
Leg length	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Hip circumference guidelines	32.5"–44.5" 83–112.5cm		34.5"–50" 87–127cm		38.5"–59.5" 98–150.5cm	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Maternity hip circumference guidelines: Small: ≤56" (142cm), Medium: ≤62" (157.5cm), Large: ≤70" (178cm)

Styles	Colors						
Closed Toe = ● Open Toe = ○	Midnight Blue 09	Chai 35	Pecan 40	Mulberry 60*	Espresso 89*	Graphite 91	Black 99
							
Calf (C)	●	●○	●	●	●	●	●○
Calf w/grip-top (C/S)*							●
Thigh w/grip-top (N) Sensinnov band grip-top	●	●○	●	●	●	●	●○
Pantyhose (P)	●	●○	●	●	●	●	●○
Maternity Pantyhose (M)**							●

* 20–30mmHg ONLY

** Maternity styles are only available in 15–20mmHg and 20–30mmHg. They are made to order and shipped in 7–10 business days

Product code example: 843NMLW09					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
84	3	N	ML	W	09

Product code example: 841MLSW99					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
84	1	M	LS	W	99

Motion

CUSHIONED COTTON 360

Motion Cushioned Cotton is soft and breathable with an extra-padded sole for optimal comfort, and is ideal for sensitive skin.

Compression level

20–30mmHg

Features and benefits

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Flat, low-friction toe seam

Fiber content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Work shoes and boots
- Casual shoes and sports shoes



362CMLM99
20–30mmHg, Closed Toe Calf
Medium Long in Black





362CMLM00
20–30mmHg, Closed Toe Calf
Medium Long in White

Motion Cushioned Cotton 360 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm		11.5"–14" 29–36cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" 35.5–51cm		16"–24" 41–61cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

Style	Colors	
Closed Toe = ●	White 00	Black 99
		
Calf (C)	●	●

Product code example: 362CLSW00					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
36	2	C	LS	W	00

Motion

HIGH TECH 412

Motion High Tech delivers optimal comfort where it counts with an anatomically-contoured shape and a cushioned sole.

Compression level

20–30mmHg

Features and benefits

- High-tech fibers for thermal control
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Additional compression in calf can improve circulation

Fiber content

- 64% Nylon, 17% Polyester, 19% Spandex
- Non-latex

Best suited for

- Management of early stages of Chronic Venous Disease (CVD)
- People who stand on their feet all day
- Walking and everyday activities
- Casual and occupational footwear



412CML59
20–30mmHg, Closed Toe Calf
Medium Large in Red



412CMM62
20–30mmHg, Closed Toe Calf
Medium Medium in Steel Blue

Motion High Tech 412 sizing chart

	Small			Medium			Large			Extra Large		
Ankle circumference	7"-8.5" 18-21.5cm			8.5"-10" 21.5-25.5cm			10"-11.5" 25.5-29cm			11.5"-13" 29-33cm		
Calf circumference	11"-15.5" 28-39cm			12.5"-17.5" 31.5-45cm			14"-20" 35.5-51cm			16"-23" 41-58cm		
Women's shoe size	4.5-7.5	8-10.5	11-13.5	4.5-7.5	8-10.5	11-13.5	8-10.5	11-13.5	>14	8-10.5	11-13.5	>14
Men's shoe size	≤6.5	7-9.5	10-12.5	≤6.5	7-9.5	10-12.5	7-9.5	10-12.5	13-15.5	7-9.5	10-12.5	13-15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

Style		Colors						
Closed Toe = ●		White 00	Blue 50	Lime 54	Pink 56	Red 59	Steel Blue 62	Black 99
								
Calf (C)		●	●	●	●	●	●	●

Product code example: 412CML00				
Series	mmHg Level	Style	Size	Color
41	2	C	ML	00

Motion

THERMOREGULATING WOOL 422

Motion Thermoregulating Wool combines the benefits of natural Merino wool with cushioned soles, making it ideal for wherever the road takes you!

Compression level

20–30mmHg (also available in 15–20mmHg, see p. 69)

Features and benefits

- Merino wool is softer and thinner than regular wool
- Wool fibers naturally help regulate body temperature and transport sweat away from the skin
- Natural properties of wool make these socks odor resistant
- Cushioned, padded soles may prevent blisters

Fiber content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Non-latex

Best suited for

- Management of symptoms related to Chronic Venous Disease (CVD)
- People who stand on their feet all day
- Work boots, steel-toed shoes, and other occupational footwear
- Adventurers, hikers, campers, craftsman, and gardeners



422CML31
20–30mmHg, Closed Toe Calf
Medium Large in Olive



422CLL12
20–30mmHg, Closed Toe Calf
Large Large in Charcoal

Motion Thermoregulating Wool 422 sizing chart

	Small			Medium			Large		
Ankle Circumference	7"–8.5" 18–21.5cm			8.5"–10" 21.5–25.5cm			10"–11.5" 25.5–29cm		
Calf Circumference	11"–15.5" 28–39cm			12.5"–17.5" 31.5–45cm			14"–20" 35.5–51cm		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

Style

Colors

Closed Toe = ●

Charcoal
12

Olive
31

Calf (C)

●

●

Product code example: 422CSL12				
Series	mmHg Level	Style	Size	Color
42	2	C	SL	12

sigvaris

Specialty

NATURAL RUBBER 500

Natural Rubber offers patients a durable, long-term containment solution for compression therapy with minimal fiber fatigue throughout the day.

Compression levels

30–40mmHg • 40–50mmHg
50–60mmHg

Features and benefits

- Ideal for patients with medical concerns
- Thigh band with low profile grip-top
- Extended sizing for a custom-like fit

Fiber content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- Caution: contains natural rubber latex that may cause allergic reaction

Best suited for

- Varicosities
- Chronic edema
- Venous stasis
- Early stages of lymphedema
- Active venous ulcers



504CM3077
40–50mmHg, Open Toe Calf
Medium Full Short in Beige



Silicone dotted grip-top band
available for thigh

Specialty Natural Rubber 500 sizing chart

	Small				Medium				Large				Extra Large			
Ankle circumference	7.5"–9" 19–23cm				9"–12" 23–30cm				12"–15.5" 30–40cm				15.5"–18" 40–46cm			
Calf circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	11"–15" 28–38cm		<17" <43cm		12.5"–17" 32–43cm		<19" <48cm		14.5"–19" 36–48cm		<22.5" <57cm		18"–22.5" 45–57cm		<27" <68cm	
Calf circumference for 505C (Full only)	505C: full calf only		12.5"–17" 32–43cm		505C: full calf only		14.5"–19" 36–48cm		505C: full calf only		18"–22.5" 45–57cm		505C: not available in XL			
Calf length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm	≤15.5" 40cm	>15.5" 40cm
Thigh circumference	Average		Full		Average		Full		Average		Full		Custom sizes available in Natural Rubber, cotton, and opaque styles. Please contact Customer Care for details.			
	16"–21" 41–54cm		<24" 61cm		19"–24" 48–61cm		<28" 71cm		23"–28" 58–71cm		<32" 81cm					
Thigh circumference for 505T (Full only)	505T: full only		<24" 61cm		505T: full only		<28" 71cm		505T: full only		<32" 81cm					
Leg length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	<28" 71cm	>28" 71cm	<28" 71cm	>28" 71cm	<28" 71cm	>28" 71cm	<28" 71cm	>28" 71cm	<28" 71cm	>28" 71cm	<28" 71cm	>28" 71cm				
Hip circumference guidelines	<48" 122cm				<54" 137cm				<61" 155cm							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

Styles

Color

Open Toe = **O**

Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 75 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.

Beige
77



Calf (C)	Calf only available in S, M, and L full circumference for 50–60mmHg	O
Thigh w/out grip-top (T)		O
Thigh w/grip-top (N) Silicone band grip-top	Not available in 50–60mmHg	O
Pantyhose (P)	Pantyhose only available in 30–40mmHg and only in full thigh circumference	O
Thigh w/Waist Attachment (W) Indicate leg LEFT (/L) or RIGHT (/R)	Thigh w/ waist attachment is available in average and full thigh circumferences in 30–40mmHg and in full thigh circumferences only in 40–50mmHg	O

Product code example: 504CS2O77

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
50	4	C	S2	O	77

Product code example: 503WL1O77/L

Series	mmHg Level	Style	Size	Gender / Open Toe	Color	Leg
50	3	W	L1	O	77	/L

Specialty

SECURE 550

Secure combines compression and containment with comfort and durability that's optimal for comfortably managing advanced venous edema, lymphedema, and wound care.

Compression level

20–30mmHg • 30–40mmHg
40–50mmHg

Features and benefits

- More graduated compression and containment than other traditional garments
- Durable, double-covered inlay yarns allow comfortable all-day wear

Fiber content

- 78% Nylon, 22% Spandex (Calf)
- 76% Nylon, 24% Spandex (Thigh)
- Non-latex

Best suited for

- Chronic Venous Insufficiency (CVI)
- Lymphedema that is well-managed
- Lipedema
- Varicosities
- Deep vein thrombosis (DVT)
- Superficial phlebitis
- Post-thrombotic syndrome
- Orthostatic hypotension
- People with sensitivity to latex



554CL3W99/S
40–50mmHg, Closed Toe Calf w/grip-top
Large Full Short in Black



Silicone dotted grip-top band
available for calf and thigh

Specialty Secure 550 Sizing Chart

	Medium				Large				X Large				XX Large			
Ankle circumference	8.5"-10" 21.5-25.5cm				10"-11.5" 25.5-29cm				11.5"-13" 29-33cm				13"-14.5" 33-37cm			
Calf circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	12"-17.5" 31-45cm		17.5"-23" 45-59cm		13.5"-19" 35-49cm		19"-25" 49-63cm		16"-22" 41-56cm		22"-28" 56-71cm		17.5"-23" 45-59cm		23"-29" 59-74cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	20"-27" 51-69cm		23.5"-31" 60-78cm		24"-30" 61-77cm		26.5"-34" 68-86cm		27"-33.5" 69-85cm		30"-37" 76-94cm		29.5"-36.5" 75-93cm		33"-40" 84-102cm	
Leg length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Size	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4	G1	G2	G3	G4

Women's Shoe Sizes: Short Leg 5.5-9.5, Long Leg 8-12 Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

Styles

Colors

Closed Toe = ●
Open Toe = ○

Beige
77

Black
99

	Beige 77		Black 99	
Calf w/grip-top (C/S) Silicone band grip-top	●○		●○	
Thigh w/grip-top (N) Silicone band grip-top	●○		●○	

Product code example: 553NX2W77

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
55	3	N	X2	W	77

Specialty

SECURE ARMSLEEVE, GLOVE, AND GAUNTLET 560

The Secure Armsleeve combines compression and containment with comfort and durability that's optimal for managing lymphedema.

Compression level

15–20mmHg • 20–30mmHg

Features and benefits

- More graduated compression and containment than traditional garments
- Full-fingered Glove and Gauntlet available
- Provides UPF 50+ sun protection

Fiber content

- 78% Nylon, 22% Spandex (Armsleeve)
- 60% Nylon, 40% Spandex (Glove, Gauntlet)
- Non-latex

Best suited for

- Chronic edema
- Lymphedema that is well-managed
- Post-operative edema
- Post-mastectomy care



**Silicone dotted
grip-top band**



562GM78
Gauntlet Medium in Beige



562FM78
Glove Medium in Beige

Specialty Secure Armsleeve 560 sizing chart

	Small				Medium				Large				X Large				XX Large			
Wrist	5.5"–6.5" 14–16cm				6.5"–7" 16–18cm				7"–8" 18–20cm				8"–9" 20–23cm				9"–10" 23–26cm			
Elbow	Average		Full		Average		Full		Average		Full		Average		Full		Average		Full	
	8"–10" 20–26cm		10"–13.5" 26–34cm		9"–11.5" 23–29cm		11.5"–14.5" 29–37cm		10"–12.5" 26–32cm		12.5"–15.5" 32–40cm		11.5"–13.5" 29–34cm		13.5"–16.5" 34–42cm		12.5"–14.5" 32–37cm		14.5"–17.5" 37–45cm	
Upper arm	Average		Full		Average		Full		Average		Full		Average		Full		Average		Full	
	10"–13.5" 26–34cm		13.5"–15.5" 34–40cm		11.5"–14.5" 29–37cm		14.5"–17" 37–43cm		12.5"–15.5" 32–40cm		15.5"–18" 40–46cm		13.5"–16.5" 34–42cm		16.5"–19" 42–48cm		14.5"–17.5" 37–45cm		17.5"–20.5" 45–52cm	
Arm length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4	G1	G2	G3	G4

Specialty Secure Glove and Gauntlet 560 sizing chart

	X Small		Small		Medium		Large		X Large		XX Large	
Palm circumference	6"–6.5" 16–17cm		7"–7.5" 18–19cm		8"–8.5" 20–21cm		8.5"–9" 22–23cm		9.5"–10" 24–25cm		10.5"–11" 26–28cm	
Size	XS		S		M		L		XL		G	

Styles

Colors

Sleeve, Glove, /
and Gauntlet = ○

Blue
50


Dusty
Rose
65

Beige
78

Cocoa
84

Black
99

	Blue 50	Dusty Rose 65	Beige 78	Cocoa 84	Black 99
Armsleeve (A)		○	○	○	○
Glove (F)	○	○	○	○	○
Gauntlet (G)	○	○	○	○	○



Product code example: 561AL278

Series	mmHg Level	Style	Size	Color
56	1	A	L2	78

Product code example: 562FM78

Series	mmHg Level	Style	Size	Color
56	2	F	M	78

Product code example: 562GM78

Series	mmHg Level	Style	Size	Color
56	2	G	M	78

sigvaris

Specialty

SECURE LITE

ARMSLEEVE 570

The Secure Lite Armsleeve combines ultra-soft materials with a seamless design to manage and prevent lymphedema with maximum comfort.

Compression level

15–20mmHg • 20–30mmHg

Features and benefits

- Soft, comfortable materials that glide on and off with ease
- Micro-dot grip-top band to ensure garment stays in place all day long
- Provides UPF 50+ sun protection

Fiber content

- 69% Nylon, 31% Spandex (Armsleeve)
- 60% Nylon, 40% Spandex (Glove, Gauntlet)
- Non-latex

Best suited for

- Chronic edema
- Lymphedema that is well-managed
- Post-operative edema
- Post-mastectomy care



Micro-dot grip-top band



562GM65

Gauntlet Medium in Dusty Rose



562FM65

Glove Medium in Dusty Rose

Specialty Secure Lite Armsleeve 570 Sizing Chart

	Small	Medium	Large	X Large	XX Large
Wrist	5.5"–6.5" 14–16cm	6.5"–7" 16–18cm	7"–8" 18–20cm	8"–9" 20–23cm	9"–10" 23–26cm
Elbow	8"–13.5" 20–34cm	9"–14.5" 23–37cm	10"–15.5" 26–40cm	11.5"–16.5" 29–42cm	12.5"–17.5" 32–45cm
Upper Arm	10"–15.5" 26–40cm	11.5"–17" 29–43cm	12.5"–18" 32–46cm	13.5"–19" 34–48cm	14.5"–20.5" 37–52cm
Size	S	M	L	X	G

Designed to fit arm length range of 14"–23" (36cm–59cm)

Specialty Secure Glove and Gauntlet 560 Sizing Chart


	X Small	Small	Medium	Large	X Large	XX Large
Palm Circumference	6"–6.5" 16–17cm	7"–7.5" 18–19cm	8"–8.5" 20–21cm	8.5"–9" 22–23cm	9.5"–10" 24–25cm	10.5"–11" 26–28cm
Size	XS	S	M	L	XL	G

Styles Colors

Sleeve, Glove, and Gauntlet = O	Dusty Rose Floral 16*	Blue 50**	Dusty Rose 65*	Blue Chevron 70**	Beige 78	Cocoa 84	Black 99
Armsleeve (A)	O			O	O	O	O
Glove (F)		O	O	O	O	O	O
Gauntlet (G)		O	O	O	O	O	O

* Dusty Rose Glove and Gauntlet match the Dusty Rose Floral Armsleeve

** Blue Glove and Gauntlet match the Blue Chevron Armsleeve



Product code example: 571AX16

Series	mmHg Level	Style	Size	Color
57	1	A	X	16

Product code example: 562FM70

Series	mmHg Level	Style	Size	Color
56	2	F	M	70

Product code example: 562GM70

Series	mmHg Level	Style	Size	Color
56	2	G	M	70

Specialty

ADVANCE ARMSLEEVE 910

The Advance Armsleeve combines ultimate softness with silver ion properties and a Sensinnov grip-top band that's ideal for all-day comfort and effectiveness.

Compression level

20–30mmHg

Features and benefits

- Ideal for patients with arm lymphedema or post-mastectomy care
- Sensinnov top band ensures armsleeve stays in place all day long
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable to provide comfortable, year-round wear

Fiber content

- 50% Nylon, 28% Rayon, 22% Spandex
- Non-latex

Best suited for

- Chronic edema
- Early stages of lymphedema



Sensinnov grip-top band

Specialty Advance Armsleeve 910 Sizing Chart


	Small		Medium				Large				X Large			
Palm	6.5"–8" 17–21cm		7"–8.5" 18–22cm				7.5"–9" 19–23cm				8"–9.5" 20–24cm			
Wrist	5.5"–6.5" 14–16cm		6.5"–7" 16–18cm				7"–8" 18–20cm				8"–9" 20–23cm			
Upper Mid-Arm	Regular		Regular		Plus		Regular		Plus		Regular		Plus	
	10"–13.5" 26–34cm		11.5"–14.5" 29–37cm		14.5"–17" 37–43cm		12.5"–15.5" 32–40cm		15.5"–18" 40–46cm		13.5"–16.5" 34–42cm		16.5"–19" 42–48cm	
Arm Length	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm	≤17" 43cm	>17" 43cm
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL	LR/P	LL/P	XR	XL	XR/P	XL/P

Style

Sleeve = ●

Beige
77

Armsleeve w/gauntlet (G)
Sensinnov band grip-top



Armsleeve w/ Gauntlet
product code example: 912ALRG77/S

Series	mmHg Level	Style	Size	w/ Gauntlet	Color	Grip-top
91	2	A	LR	G	77	S

Armsleeve w/ Gauntlet
product code example: 912AMLG77/PS

Series	mmHg Level	Style	Size	w/ Gauntlet	Color	Plus size	Grip-top
91	2	A	ML	G	77	/P	S

Specialty

DIABETIC COMPRESSION SOCKS 602

Diabetic Compression Socks are the only patented, clinically-proven graduated compression socks constructed for people with diabetes.

Compression level

18–25mmHg (also available in 8–15mmHg, see p. 70)

Features and benefits

- Padded soles and invisible toe seam add comfort and protection
- Compression-free foot for optimum comfort
- Double-covered inlay yarns provide comfort and allow product to easily glide on and off
- *drirelease*® yarns add moisture wicking properties

Fiber content

- 76% *drirelease* (65% Polyester, 11% Cotton), 19% Nylon, 5% Spandex
- Non-latex

Best suited for

- Individuals with diabetes and symptoms related to Chronic Venous Disease (CVD)
- Individuals with peripheral neuropathy
- Patients with an ABI of >0.6



602CMLM00
18–25mmHg, Closed Toe Calf
Medium Long in White



Padded soles and invisible toe seam
add comfort and protection

Specialty Diabetic Compression Socks 602 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"–8.5" 18–21.5cm		8.5"–10" 21.5–25.5cm		10"–11.5" 25.5–29cm		11.5"–14" 29–36cm	
Calf circumference	11"–15.5" 28–39cm		12.5"–17.5" 31.5–45cm		14"–20" (5.5–51cm)		16"–24" 41–61cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL	XS*	XL*

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14
*XS & XL Only Available In Men's Sizing

Style	Color
<div>Closed Toe = ●</div>	<div>White 00</div> <div></div>
<div>Calf (C)</div>	<div>●</div>

Product code example: 602CSLW00					
Series	mmHg Level	Style	Size	Gender / Open Toe	Color
60	2	C	SL	W	00

Measuring Guide for Lower Inelastic Compression

Patient Name: _____ Phone #: _____ Date: _____

Patient Address: _____ City: _____ State: _____ Zip: _____

Prescribing Doctor's Name (Printed): _____

NPI #: _____ Phone #: _____

Doctor Signature: _____ ICD 10 Diagnosis Code: _____

Body Part*
☐ Calf ☐ Knee ☐ Foot

Inelastic Wraps
☐ Compreflex ☐ Coolflex ☐ Compreboot
☐ Light ☐ Standard ☐ Plus

Product style types may vary, please see the 2024 Sigvaris Product Catalog for Inelastic Wraps, Lymphedema Wear, Bandages and Therapy Supplies for more information or call the US Customer Care Solution Center. Available Monday–Friday, 8am–7pm EST at 1.800.322.7744, ext. 1.

Inelastic Wrap: Lower Extremity Measurements

Circumference
Left Right

F _____

E _____

D _____

C1 _____

C _____

B _____

A _____

I _____

J _____

H Length

H _____

K _____

Gluteal Fold

Top of Patella

Fibular Head

Lateral Malleolus
Outer ankle bone

Ø Point

K: Measure **medial** length from heel to 1st metatarsal head

Inelastic Wrap: Compreflex Calf Garment Measurements

Circumference
Left Right

C _____

A _____

G Length
Regular: 30–33.9cm
Tall: 34–42cm

G _____

Fibular Head

Largest Calf

G Length
Lateral malleolus to fibular head

Largest Ankle
Lateral Malleolus
Outer ankle bone

K: Measure **medial** length from heel to 1st metatarsal head

Note: on the following pages, these icons indicate options for time of wear.

Indicates Day wear

Indicates Night wear

For a PDF of the product ordering form, scan the QR code or visit sigvaris.group/mceforms

*Measure point letter may vary with inelastic garments, please see the 2024 Sigvaris Product Catalog for Inelastic Wraps, Lymphedema Wear, Bandages and Therapy Supplies for a product specific point.



Precise measurements are essential to guarantee the indicated level of compression and to ensure a proper fit. The measuring instructions included in the following pages are designed for your understanding.

The measurement chart and graphic for each product is also included on corresponding product specific order form that can be found online at: sigvaris.group/en-us/expertise/download-center

For additional guidance with measuring or fitting of a garment, please contact the Customer Care Solution Center, available Monday–Friday, 8am–7pm at **1.800.322.7744**, ext. 1.

CALF & FOOT

Circumferences:

For Compreflex Calf garments, measure circumferences where indicated on the measurement chart at the largest ankle and largest calf. Measurements may be recorded on the corresponding lines A and C.

For all other Calf garments, measure circumferences where indicated on the measurement chart, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone), continue to the 15cm, 25cm, and 30cm dots. Measurements may be recorded on the corresponding lines A, B, C, C (left/right).

Length:

Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Fibular Head**. Measurements may be recorded on the box G.

Foot:

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.



Measurement tips

- 1 When working with a patient tasked with self-measuring, it is best to provide the measuring graphic (included in this guide) as a visual guide. Studies also show that the use of medical terminology may be confusing and even frightening for the patient. For the best result, limit the use of medical terminology.
- 2 Remind your patient, the tape or measurement tool should be held taut against the skin but should not dent the skin.
- 3 If patient has a caregiver or companion, it may be best to involve them for the placement of the tape or the reading of the measurements.
- 4 Using a body marking pen or makeup pencil may also help with recording landmarks on the limb for future reference.
- 5 Do not mark on fragile or broken skin.

Compreflex Standard CALF & FOOT

20-50mmHg



Colors



Black

Beige

When to Use

- Adjustable ranges of resting compression are desired (20-50mmHg)
- Detached foot piece is desired
- No compression in the liner is desired
- Patient requires a calf and foot garment

Ideal for

- Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post-surgical edema




Compreboot® Standard

Variable compression
(see page 54 for details)



Basic Liner

Non-Compressive
(see page 53 for details)



Accutab® Devices

Accutab Devices
(visit website for details)

Compreflex Standard CALF

20-50mmHg



Colors



Black

Beige

When to Use

- Varying ranges of adjustable resting compression are desired (20-50mmHg)
- No compression in the liner is desired
- No Compreboot needed
- Patient requires compression in the calf only or a different foot option (size or boot)
- Ability to wear standard footwear over liner is desired


Ideal for

- Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post-surgical edema



Basic Liner

Non-Compressive
(see page 53 for details)



Accutab® Devices

Accutab Devices
(visit website for details)

Medically Complex Edema

COMPREFLEX WRAPS




Variable compression Compreboot for added foot and ankle compression (see page 56 for details)

Compreflex Standard Calf & Foot Sizing Chart					
	Small	Medium	Large	X Large	XX Large
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm
Regular	1401-BKR	1402-BKR	1403-BKR	1404-BKR	1405-BKR
Tall	1401-BKT	1402-BKT	1403-BKT	1404-BKT	1405-BKT
Regular	1411-BKR	1412-BKR	1413-BKR	1414-BKR	1415-BKR
Tall	1411-BKT	1412-BKT	1413-BKT	1414-BKT	1415-BKT

Tall garments ship with Regular length Compreboot. If Long Compreboot is required, purchase Compreflex Standard Calf and order desired Compreboot separately.

REQUIRED MEASUREMENTS



Compreflex Standard – Calf Sizing Chart (A–C)					
	Small	Medium	Large	X Large	XX Large
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm
Regular	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
Tall	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT
Regular	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
Tall	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT



For a PDF of the product ordering form, scan the QR code or visit sigvaris.group/mceforms


Compreflex Transition

CALF

20-50mmHg



Colors



When to Use

- Patient requires a compressive liner with 10-15mmHg in the foot and ankle
- Ability to wear standard footwear over liner is desired

Ideal for

- Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post-surgical edema



Transition
Liner

10-15 mmHg of mild compression in the foot and ankle
(see page 53 for details)



Accutab Devices

Accutab Devices
(visit website for details)


Compreflex Complete

CALF

20-50mmHg



Colors



When to Use

- Patient requires a compressive liner with 10-15mmHg in the foot, ankle, and calf
- Ability to wear standard footwear over liner is desired


Ideal for

- Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post-surgical edema



Complete
Liner

10-15 mmHg of mild compression in foot, ankle, and calf
(see page 53 for details)




Accutab Devices

Accutab Devices
(visit website for details)

Medically Complex Edema

COMPREFLEX WRAPS

REQUIRED MEASUREMENTS



Compreflex Complete Calf
20-50mmHg
Medium Regular in Beige
1412-FC-BKR



For a PDF of the product ordering form, scan the QR code or visit sigvaris.group/mceforms

Compreflex Transition – Calf Sizing Chart (A–C)

		Small	Medium	Large	X Large	XX Large
	C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
	A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm
Black	Regular	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
	Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT
Beige	Regular	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
	Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT

Compreflex Complete – Calf Sizing Chart (A–C)

		Small	Medium	Large	X Large	XX Large
	C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
	A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm
Black	Regular	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
	Tall	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT
Beige	Regular	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
	Tall	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT

Compreflex Strap Extender

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB

- Provides up to 10cm of additional circumferential size range to Compreflex calf wraps (excludes Compreflex Contain)
- Detaches easily
- For patients outside of standard size ranges of Compreflex calf wraps



Compreflex Contain

CALF

20–50mmHg



Colors



Black



Beige

When to Use

- Patient requires a compressive liner with 10–15mmHg in the foot and ankle
- Ability to wear standard footwear over liner is desired


Ideal for

- Lymphedema
- Edema
- Lipedema

Features and benefits

- Thicker material provides more containment than other Compreflex products to address lower extremity lymphedema
- Garment features a moderately-contoured calf design to fit edematous legs
- 20–50mmHg of adjustable, gradient, inelastic compression
- Stretch panel holds garment in place for easy application and high compliance
- Transition Liners feature an antimicrobial wash for bacteriostatic properties and odor control and an average of 10–15mmHg of compression in the foot and ankle (see page 53 for details)

Billable under
Medicare code **A6583**



Transition Liner

10–15 mmHg of mild
compression in the foot
and ankle
(see page 53 for details)

Fiber content

- 48% Cotton • 25% Polyamide
13% Polyester • 9% Polyurethane Foam
5% Elastane
- Non-latex

Best suited for

- Lymphedema
- Edema
- Lipedema



Compreflex Contain (bottom) is thicker
than other Compreflex products (top)




Also available in Black

Medically Complex Edema

COMPREFLEX WRAPS

REQUIRED MEASUREMENTS



G Length

Regular:
30–33.9cm
Tall:
34–42cm

Compreflex Contain – Calf Sizing Chart (A–C)

	Small	Medium	Large	X Large	XX Large
C	29–39cm	34–44cm	40–50cm	46–56cm	54–64cm
A	19–26cm	21–30cm	26–36cm	31–41cm	36–46cm
Black					
Regular	3401-UC-BKR	3402-UC-BKR	3403-UC-BKR	3404-UC-BKR	3405-UC-BKR
Tall	3401-UC-BKT	3402-UC-BKT	3403-UC-BKT	3404-UC-BKT	3405-UC-BKT
Beige					
Regular	3411-UC-BKR	3412-UC-BKR	3413-UC-BKR	3414-UC-BKR	3415-UC-BKR
Tall	3411-UC-BKT	3412-UC-BKT	3413-UC-BKT	3414-UC-BKT	3415-UC-BKT

Identifying the optimal product size ensures a proper fit and the correct distribution of compression.

Compreflex Contain Strap Extender

- Provides up to 10cm of additional length to a single strap
- 1 unit extends 1 strap
- Detaches easily
- For patients outside of standard size range of Compreflex Contain calf wrap
- This Strap Extender works only with the Compreflex Contain wrap

	1 Unit	10 Unit Pack
Black	3401-SE1	3401-SEB
Beige	3411-SE1	3411-SEB



For a PDF of the product ordering form, scan the QR code or visit sigvaris.group/mceforms

Coolflex
Standard

CALF

20–50mmHg



Color



Black

When to Use

- Varying ranges of adjustable resting compression are desired (20–50mmHg)
- No compression in the liner is desired
- Simplified one-handed or lateral donning is desired
- More breathable fabric is desired
- Desire to select separate foot piece

Ideal for

- Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Fragile, sensitive skin
- Well-managed Lymphedema
- Lipedema

One-Handed Donning

Spacer Fabric is highly breathable for cool-wearing comfort


Durable metal hook closures for maximum strength and durability

Basic Liner

Non-Compressive
(see page 53 for details)



Coolflex Standard Calf Sizing Chart (A–C1)					
SIZE	Small	Small Max	Medium/Large	Medium / Large Max	X Large / XX Large
C1	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm
C	28–39cm	35–50cm	34–50cm	45–64cm	45–64cm
B	23–33cm	27–38cm	29–43cm	35–49cm	37–54cm
A	16–26cm	16–26cm	23–33cm	23–33cm	30–42cm
Left					
Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
Right					
Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R



Basic Liner

Non-compressive

Features

- Non-compressive liner for the protection of the skin
- Antimicrobial wash for bacteriostatic properties and odor control
- Flat seams for added comfort
- Available in Heathered Silver color only

Size: One size fits all	
Pair	10 Pair Pack
1700 - NK1	1700 - NKC



- 13% Nylon
- 1% Spandex
- 13% Cotton
- 73% Polyester

Transition Liner

10–15mmHg

Features

- Mild compression in foot and ankle only
- Antimicrobial wash for bacteriostatic properties and odor control
- Extra thick sole from toe to heel for added comfort
- Available in Black and Beige

Transition Liner Sizing Chart						
SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19–26cm	29–39cm	<7½	<6½	1601-CK	1611-CK
Medium	21–30cm	34–44cm	7.5–9.5	6.5–8.5	1602-CK	1612-CK
Large	26–36cm	40–50cm	10–12	9–11	1603-CK	1613-CK
X Large	31–41cm	46–56cm	12.5–15	11.5–14	1604-CK	1614-CK
XX Large	36–46cm	54–64cm	12.5–15	11.5–14	1605-CK	1615-CK



- 16% Nylon
- 4% Spandex
- 12% Cotton
- 68% Polyester

Complete Liner

10–15mmHg

Features

- Mild compression in foot and ankle with graduated compression in calf
- Antimicrobial wash for bacteriostatic properties and odor control
- Extra thick sole from toe to heel for added comfort
- Available in Black and Beige

Complete Liner Sizing Chart						
SIZE	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe	Pair (Black)	Pair (Beige)
Small	19–26cm	29–39cm	<7½	<6½	1601-CK-F1	1611-CK-F1
Medium	21–30cm	34–44cm	7.5–9.5	6.5–8.5	1602-CK-F1	1612-CK-F1
Large	26–36cm	40–50cm	10–12	9–11	1603-CK-F1	1613-CK-F1
X Large	31–41cm	46–56cm	12.5–15	11.5–14	1604-CK-F1	1614-CK-F1
XX Large	36–46cm	54–64cm	12.5–15	11.5–14	1605-CK-F1	1615-CK-F1



- 19% Nylon
- 7% Spandex
- 11% Cotton
- 63% Polyester

Compreboot Lite Foot

Variable compression



Features

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild therapeutic compression

Color



Black



Compreboot Standard Foot

Variable compression



Features

- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape

Colors



Black



Beige



Compreboot Plus Foot

Variable compression



Features

- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate therapeutic compression

Color



Black



Coolflex Standard Foot

Variable compression



Features

- Low profile design allows boot to fit in most footwear*
- Spacer Fabric is highly breathable for cool-wearing comfort
- Durable metal clasps for closure

Color



Black



Foot compression garments are pictured with the Basic Liner. Basic Liner is not included.

* Data on file.

Medically Complex Edema

FOOT COMPRESSION WEAR

Circumference


Left Right

I _____

J _____

Length

K _____



K: Measure **medial** length from heel to 1st metatarsal head

Measuring Instructions

Foot Size:

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure length from heel to 1st metatarsal head

Compreboot Lite Foot Size Chart (I-K)
(Available in Black Only)

	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤33cm	≤33cm
K	14-18cm	19-23cm	20-22cm	25-27cm
Black	1803-BLR	1803-BLL	1805-BLR	1805-BLL

Compreboot Standard Foot Size Chart (I-K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL

Compreboot Plus Foot Size Chart (I-K)
(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-34cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-25cm	≤30cm	≤28cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL

Coolflex Standard Foot Size Chart (I-K)
(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-38cm	28-38cm	30-43cm	30-43cm	34-49cm	34-49cm
J	22-26cm	22-26cm	25-30cm	25-30cm	28-33cm	28-33cm
K	14-18cm	19-21cm	18-20cm	21-23cm	20-22cm	23-25cm

Left

Black	2101-BTR-L	2101-BTL-L	2103-BTR-L	2103-BTL-L	2105-BTR-L	2105-BTL-L
--------------	------------	------------	------------	------------	------------	------------

Right

Black	2101-BTR-R	2101-BTL-R	2103-BTR-R	2103-BTL-R	2105-BTR-R	2105-BTL-R
--------------	------------	------------	------------	------------	------------	------------

Compreflex Strap Extender

- Provides up to 10cm of additional circumferential size range to Compreflex calf wraps (excludes Compreflex Contain)
- Detaches easily
- For patients outside of standard size ranges of Compreflex calf wraps

	1 Unit	10 Unit Pack
Black	1401-SE1	1401-SEB
Beige	1411-SE1	1411-SEB



For a PDF of the product ordering form, scan the QR code or visit sigvaris.group/mceforms

970 Opaque WOMEN & MEN

Dynaven socks and hosiery provide opaque coverage and support for legs following a vein procedure or for everyday wear.

Compression levels

15–20mmHg • 20–30mmHg
30–40mmHg

Fiber content

- 75% Nylon, 25% Spandex
(content may vary slightly by style)
- Not made with natural rubber latex

Features

- Glide on & off easily
- Reinforced heel area
- Roomy toe area
- Wide top bands



**Silicone dotted
grip-top band**
available for calf
and thigh

920 Opaque MEN

Dynaven opaque ribbed socks provide support for legs following a vein procedure or for everyday wear.

Compression levels

15–20mmHg • 20–30mmHg
30–40mmHg

Fiber content

- 80% Nylon, 20% Spandex
- Not made with natural rubber latex

Features

- Ribbed sock style
- Glide on & off easily
- Roomy toe area
- Wide top bands



**Silicone dotted
grip-top band**
available for calf

Dynaven Opaque 970 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"-8.5" 18-21.5cm		8.5"-10" 21.5-25.5cm		10"-11.5" 25.5-29cm		11.5"-13" 29-33cm	
Calf circumference	11"-15.5" 28-39cm		12.5"-17.5" 31.5-45cm		14"-20" 35.5-51cm		16"-23" 40.5-58cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"-26" 45.5-66cm		19"-28.5" 48-73cm		20.5"-33" 52-84cm		24"-36" 61-92cm	
Leg length	Short	Long	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Hip circumference	<45.5" <115cm		<51" 129.5cm		<60.5" 154cm		N/A	
Size	SS	SL	MS	ML	LS	LL	XS*	XL*

Women's shoe sizes: Short 5.5-9.5 Long 8-12 Men's shoe sizes: Short 7-11 Long 9.5-14

* XS & XL not available in Dynaven pantyhose styles

Styles

Colors

Closed Toe = ●
Open Toe = ○

Light
Beige
66

Black
99

970 Women's opaque
calf, calf w/grip-top,
thigh w/grip-top,
and pantyhose
(C, C/S, N, & P)

970 Men's opaque
thigh w/grip top (N)

970 Unisex opaque calf,
calf w/grip-top, thigh w/
grip-top, and pantyhose
(CO, CO/S, NO, and PO)**

** Dynaven unisex opaque pantyhose only
available in 20-30mmHg and 30-40mmHg

Product code example: 972CSLW66

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
97	2	C	SL	W	66

Dynaven Opaque 920 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"-8.5" 18-21.5cm		8.5"-10" 21.5-25.5cm		10"-11.5" 25.5-29cm		11.5"-13" 29-33cm	
Calf circumference	11"-15.5" 28-39cm		12.5"-17.5" 31.5-45cm		14"-20" 35.5-51cm		16"-23" 40.5-58cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL	XS	XL

Men's shoe sizes: Short 7-11 Long 9.5-14

Styles

Colors

Closed Toe = ●

Light
Beige
66

Black
99

920 Men's opaque ribbed
calf and calf w/grip-top
(C & C/S)

Product code example 922CMLM99

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
92	2	C	ML	M	99

980 Sheer WOMEN

Dynaven sheer hosiery offers medical grade compression and features sheer fabric for a fashionable look. Ideal for everyday wear, Dynaven sheer hosiery supports leg and vein health.

Compression levels

15–20mmHg • 20–30mmHg

Fiber content

- 65% Nylon, 35% Spandex (content may vary slightly by style)
- Not made with natural rubber latex

Features

- Sheer fabric for a fashionable look
- Glide on & off easily
- Lace silicone grip-top band for thigh w/grip-top styles
- Cotton-lined gusset for pantyhose styles
- Wide top bands



Silicone grip-top band
available for thigh

960 Cushioned WOMEN & MEN

Dynaven Cushioned socks feature cushioned heel, sole, and toe areas. Ideal for everyday wear or following a vein procedure, the Dynaven Cushioned sock gives legs the support and comfort they need to go the distance.

Compression levels

15–20mmHg • 20–30mmHg

Fiber content

- 82% Nylon, 18% Spandex
- Not made with natural rubber latex

Features

- Glide on & off easily
- Cushioned heel, sole, & toe areas
- Wide top bands
- Single sizing system for both men and women



Dynaven Sheer 980 sizing chart

	Small		Medium		Large	
Ankle circumference	7"-8.5" 18-21.5cm		8.5"-10" 21.5-25.5cm		10"-11.5" 25.5-29cm	
Calf circumference	11"-15.5" 28-39cm		12.5"-17.5" 31.5-45cm		14"-20" 35.5-51cm	
Calf length	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Thigh circumference	18"-26" 45.5-66cm		19"-28.5" 48-73cm		20.5"-33" 52-84cm	
Leg length	Short	Long	Short	Long	Short	Long
	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm	≤28.5" 72cm	>28.5" 72cm
Hip circumference	<45.5" 115cm		<51" 129.5cm		<60.5" 154cm	
Size	SS	SL	MS	ML	LS	LL

Women's shoe sizes: Short 5.5-9.5 Long 8-12

Styles


Colors

Closed Toe = ●
Open Toe = ○

Light
Beige
77

Black
99

980 Women's Sheer calf,
thigh-high w/grip-top,
and pantyhose
(C, N, & P)*



* Dynaven sheer pantyhose only available
in closed toe

Product code example: 982NSSW77

Series	mmHg Level	Style	Size	Gender / Open Toe	Color
98	2	N	SS	W	77

Dynaven Cushioned 970 sizing chart

	Small		Medium		Large		Extra Large	
Ankle circumference	7"-8.5" 18-21.5cm		8.5"-10" 21.5-25.5cm		10"-11.5" 25.5-29cm		11.5"-13" 29-33cm	
Calf circumference	11"-15.5" 28-39cm		12.5"-17.5" 31.5-45cm		14"-20" 35.5-51cm		16"-23" 40.5-58cm	
Calf length	Short	Long	Short	Long	Short	Long	Short	Long
	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm	≤16" 41cm	>16" 41cm
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's shoe sizes: Short 5.5-9.5 Long 8-12 Men's shoe sizes: Short 7-11 Long 9.5-14

Style


Colors

Closed Toe = ●

White
00

Black
99

960 Unisex cushioned
calf (C)



Product code example 961CLS00

Series	mmHg Level	Style	Size	Color
96	1	C	LS	00

Casual Cotton

146 / 186

For women & men

Compression level

15–20mmHg

Features and benefits

- Made with premium Supima cotton for optimum softness and durability
- Flat, low-friction toe seam to prevent irritation
- Soft, stretchable, knit-in band prevents pinching

Fiber content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease



186CB00
15–20mmHg Closed Toe Calf Size B in White



146CC11
15–20mmHg Closed Toe Calf Size C in Brown



146CB30
15–20mmHg Closed Toe Calf Size B in Khaki



186CB99
15–20mmHg Closed Toe Calf Size B in Black

146 Women's Calf sizing chart

Women's shoe	5–7	7.5–9.5	10–12
Ankle circumference	7"–9" 18–23cm	7.5"–10" 19–25cm	8"–10.5" 20–27cm
Size	A	B	C

186 Men's Calf sizing chart

Men's shoe	6.5–8.5	9–11	11.5–14
Ankle circumference	7.5"–10" 19–25cm	8"–10.5" 20–27cm	8.5"–11.5" 21–29cm
Size	A	B	C

Style Colors

Closed Toe = ●

White
00Navy
10Brown
11Khaki
30Black
99

	White 00	Navy 10	Brown 11	Khaki 30	Black 99
Calf (C)	●	●	●	●	●

Microfiber Shades

143 / 183

For women & men

Compression level

15–20mmHg (also available in
20–30mmHg, see p. 22)

Features and benefits

- Advanced microfiber technology for ultimate comfort
- Glides on and off with ease
- Excellent breathability and moisture management
- Soft and durable sock for everyday wear

Fiber content

- 87% Nylon, 13% Spandex
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease



143CC02
15–20mmHg Closed Toe Calf Size C
in Mariner Stripe



183CB17
15–20mmHg Closed Toe Calf Size B
in Graphite Stripe



183CB47
15–20mmHg Closed Toe Calf Size B
in Graphite Argyle



143CA48
15–20mmHg Closed Toe Calf Size A
in Royal Blue Argyle

143 Women's Calf sizing chart

Women's shoe	5–7	7.5–9.5	10–12
Ankle circumference	7"–9" 18–23cm	7.5"–10" 19–25cm	8"–10.5" 20–27cm
Size	A	B	C

183 Men's Calf sizing chart

Men's shoe	6.5–8.5	9–11	11.5–14
Ankle circumference	7.5"–10" 19–25cm	8"–10.5" 20–27cm	8.5"–11.5" 21–29cm
Size	A	B	C

Style Colors

	Closed Toe = ●	Mariner Stripe 02	Pink Stripe 15	Graphite Stripe 17	Onyx Stripe 19	Graphite Heather 27	Black Plaid 42	Purple Argyle 46	Graphite Argyle 47	Royal Blue Argyle 48
Women's Calf (C)	●	●		●	●			●	●	●
Men's Calf (C)				●	●	●	●		●	●

Sea Island
Cotton

151 / 191

For women & men

Compression level

15–20mmHg
(also available in 20–30mmHg, see p. 12)

Features and benefits

- Knitted from the world’s rarest cotton
- Sea Island Cotton yarns are lightweight and have great strength, a silky luster, and an incredibly soft feel
- Breathable and comfortable for sensitive skin
- Ideal for travel and everyday wear

Fiber content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

151 Women’s Calf sizing chart			
Women’s shoe	5–7	7.5–9.5	10–12
Ankle circumference	7”–9” 18–23cm	7.5”–10” 19–25cm	8”–10.5” 20–27cm
Size	A	B	C

191 Men’s Calf sizing chart			
Men’s shoe	6.5–8.5	9–11	11.5–14
Ankle circumference	7.5”–10” 19–25cm	8”–10.5” 20–27cm	8.5”–11.5” 21–29cm
Size	A	B	C

Style		Colors	
Closed Toe = ●		Navy 10	Black 99
Calf (C)		●	●



All-Season Merino Wool

152 / 192

For women & men

Compression level

15–20mmHg

Features and benefits

- Soft and comfortable
- Natural properties of wool make these socks odor resistant
- Thermoregulating for year-round wear and keeps feet dry

Fiber content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease



192CB10
15–20mmHg Closed Toe Calf Size B
in Navy



192CC11
15–20mmHg Closed Toe Calf Size C
in Brown



152CC99
15–20mmHg Closed Toe Calf Size C
in Black

152 Women's Calf sizing chart

Women's shoe	5–7	7.5–9.5	10–12
Ankle circumference	7"–9" 18–23cm	7.5"–10" 19–25cm	8"–10.5" 20–27cm
Size	A	B	C

192 Men's Calf sizing chart

Men's shoe	6.5–8.5	9–11	11.5–14
Ankle circumference	7.5"–10" 19–25cm	8"–10.5" 20–27cm	8.5"–11.5" 21–29cm
Size	A	B	C

Style

Colors

Closed Toe = ●

Navy
10

Brown
11

Black
99



Calf (C)

●

●

●

Sheer Fashion

120

For women

Compression level

15–20mmHg

Features and benefits

- Fabric feels soft against the skin
- Incredible sheer design in a variety of fashionable colors
- Sheer toe and heel area
- Comfortable stretch is excellent for everyday wear
- Cotton-lined gusset on pantyhose styles

Fiber content

- 84% Nylon, 16% Spandex (content can vary by style)
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease



120PF12
15–20mmHg Closed Toe Pantyhose Size F in Charcoal



120NB29
15–20mmHg Closed Toe Thigh Size B in Taupe



120CC036
15–20mmHg Open Toe Calf Size C in Golden



120PB99
15–20mmHg Closed Toe Pantyhose Size B in Black

Styles	Colors					
Closed Toe = ● Open Toe = ○	Honey 33	Golden 36	Café 73	Taupe 29	Charcoal 12	Black 99
Calf (C)	●○	●○	●	●○	●	●
Thigh (N)	●	●	●	●	●	●
Panyhose (P)	●	●	●	●	●	●

Sheer Fashion

120 (continued)

120 Women's Calf and Thigh sizing chart			
Women's shoe	5-7	7.5-9.5	10-12
Ankle circumference	7.5"-9" 19-23cm	8"-10" 20.5-25cm	8.5"-10.5" 22-27cm
Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs., Size B: <180Lbs., Size C: <230Lbs.

Pantyhose sizing chart

[illegible]

Soft Silhouette Leggings

170

For women

Compression level

15–20mmHg

Features and benefits


- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts
- Foldable, comfort-stretch waistband provides a seamless look for all-day wearing comfort
- Advanced soft stretch microfiber fabric shapes legs and smooths skin

Fiber content


- 83% Nylon, 17% Spandex
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, family history of chronic venous disease



170LC99
15–20mmHg Legging Size C in Black



170LD89
15–20mmHg Legging Size D in Espresso



Non-binding ankle cuff
stays comfortably in place
for all day wear

170 Leggings sizing chart

Ft/Cm	Lbs/Kg	100 / 45	105 / 48	110 / 50	115 / 52	120 / 54	125 / 57	130 / 59	135 / 61	140 / 64	145 / 66	150 / 68	155 / 70	160 / 73	165 / 75	170 / 77	175 / 80	180 / 82	185 / 84	190 / 86	195 / 89	200 / 91	205 / 93	210 / 95	215 / 98	220 / 100	225 / 102	230 / 105	
4'10" / 147cm																													
4'11" / 150cm																													
5'0" / 152cm																													
5'1" / 155cm																													
5'2" / 157cm																													
5'3" / 160cm																													
5'4" / 162cm																													
5'5" / 165cm																													
5'6" / 167cm																													
5'7" / 170cm																													
5'8" / 173cm																													
5'9" / 175cm																													
5'10" / 178cm																													
5'11" / 180cm																													
6' / 183cm																													


Guidelines for Ankle Circumference

Size A: 7"–9" (18–23cm), Size B & C: 8"–10" (20–26cm), Size D: 8.5"–11.5" (22–29cm)

Style

Colors

Open Foot = ○

Midnight
Blue
09*Mulberry
60*Espresso
89*Black
99

Legging (L)

○

○

○

○

*Made on demand colors may take up to 7 days to receive

Athletic
Recovery Socks
401

For women & men

Compression level

15–20mmHg

Features and benefits

- Helps reduce lactic acid buildup¹
- Comfortable and easy to put on and take off
- *drirelease* yarns provide moisture wicking and odor control
- No compression in heel and instep for ease of donning

Fiber content

- 67% *drirelease* (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Post-workout recovery
- Mild athletic activity



401CL00
15–20mmHg Closed Toe Calf Size L
in White



401CM99
15–20mmHg Closed Toe Calf Size M
in Black

401 Calf sizing chart					
Women's shoe	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
Men's shoe	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
Ankle circumference	7"–9" 18–23cm	7.5"–10" 19–25cm	8"–10.5" 20–27cm	8.5"–11.5" 21–29cm	8.5"–11.5" 21–29cm
Size	S	M	L	X	2

Style	Colors	
Closed Toe = ●	White 00	Black 99
Calf (C)	●	●

1. Effect of Compression Stockings on Physiological Responses and Running Performance in Division III Collegiate Cross-Country Runners During a Maximal Treadmill Test. Journal of Strength and Conditioning Research. Volume 28, Number 6, June 2014. Rider, Coughlin, Hew-Butler, Goslin.

Performance
Sleeves

412V

For women & men

Compression level

20–30mmHg

Features and benefits

- Fabric resists deterioration in salt, chlorine, or fresh water
- Helps reduce lactic acid buildup¹
- Quick drying
- High-visibility reflective transfer

Fiber content

- 75% Nylon, 25% Spandex
- Non-latex

Best suited for

- Running, triathlons, endurance sports, and gym workouts
- Team sports
- Occupational wear



412VM00
20–30mmHg Sleeve Size M in White



412VM50
20–30mmHg Sleeve Size M in Blue












412VS58
20–30mmHg Sleeve Size S in Purple



412VL99
20–30mmHg Sleeve Size L in Black

412V Calf Sleeve sizing chart				
Ankle circumference	7"–8.5" 18–21.5cm	8.5"–10" 21.5–25.5cm	10"–11.5" 25.5–29cm	11.5"–13" 29–33cm
Calf circumference	11"–15.5" 28–39cm	12.5"–17.5" 31.5–45cm	14"–20" 35.5–51cm	16"–23" 41–58cm
Size	S	M	L	X

Style				
Calf Sleeve = 	White 00	Blue 50	Purple 58*	Black 99
				
Calf Sleeve (V)				

1. Effect of Compression Stockings on Physiological Responses and Running Performance in Division III Collegiate Cross-Country Runners During a Maximal Treadmill Test. Journal of Strength and Conditioning Research. Volume 28, Number 6, June 2014. Rider, Coughlin, Hew-Butler, Goslin.

* Purple not available in size X

Merino
Outdoor Sock
421

For women & men

Compression level

15–20mmHg
(also available in 20–30mmHg, see p. 30)

Features and benefits

- Soft and comfortable
- Natural properties of wool make these socks odor resistant
- Mesh flex zones for comfort and ventilation
- Cushioned, padded soles protect against blisters
- No compression in heel and instep for ease of donning
- Thermoregulating for year-round comfort

Fiber content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Non-latex

Best suited for

- Tired, heavy, restless, aching legs
- Outdoor sports and activities
- Post-workout recovery





421CS12
15–20mmHg Closed Toe Calf Size S
in Charcoal



421CX31
15–20mmHg Closed Toe Calf Size X
in Olive

421 Calf sizing chart				
Women's shoe	<7.5	7.5–9.5	10–12	12.5–15
Men's shoe	<6.5	6.5–8.5	9–11	11.5–14
Ankle circumference	7"–9" 18–23cm	7.5"–10" 19–25cm	8"–10.5" 20–27cm	8.5"–11.5" 21–29cm
Size	S	M	L	X

Style			Colors	
Closed Toe = ●			Charcoal 12	Olive 31
				
Calf (C)			●	●

Eversoft Diabetic Sock

160

For women & men

Compression level

8–15mmHg

(also available in 18–25mmHg, see p. 42)

Features and benefits

- Thick, padded soles, heels, and toes for comfort and protection
- Flat, low friction toe seam
- Soft, stretchable knit-in band prevents pinching
- Compression from ankle to knee, compression-free foot

Fiber content

- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Non-latex

Best suited for

- Individuals with diabetes
- Individuals with limited dexterity or lower tolerance for compression
- Individuals with sensitive skin



160CS00
8–15mmHg Closed Toe Calf Size S
in White



160CL00
8–15mmHg Closed Toe Calf Size L
in White

160 Calf sizing chart

Women's shoe	<7.5	7.5–9.5	10–12	12.5–15
Men's shoe	<6.5	6.5–8.5	9–11	11.5–14
Ankle circumference	7"–9" 18–23cm	7.5"–10" 19–25cm	8"–10.5" 20–27cm	8.5"–11.5" 21–29cm
Size	S	M	L	X

Style

Color

Closed Toe = ●

White
00



Calf (C)

●

ACCESSORIES

Sigvaris Accessories are created for specific medical needs, including:

- Aid in putting on and removing stockings
- Stocking hygiene to improve adherence
- Adhesives and extra liners for added comfort



ACCESSORIES

Gloves

Rubber Gloves – Ridged

The practical and protective way to apply and remove your compression stockings.



Package: 12 per Pack

Size	Product Code
Small	591R400S
Medium	591R400M
Large	591R400L
X-Large	591R400X

Package: Single Pair

Size	Product Code
Small	591RPRS
Medium	591RPRM
Large	591RPRL
X-Large	591RPRX

Glove Sizing Chart

Measure Circumference of the Hand at the Palm

6.5"–7" 16.5–18cm	7.5"–8" 19–20.5cm	8.5"–9" 21.5–23cm	9.5"–10" 24–25.5cm
Small	Medium	Large	X-Large



For more Sigvaris videos, including donning tips and product demonstrations, visit our YouTube channel: youtube.com/sigvarisgroupusa

Caution: Rubber Gloves contain natural rubber latex which may cause allergic reaction

Latex-Free Gloves

Gripping donning gloves recommended for use in a non-latex environment.



Package: 12 per Pack

Size	Product Code
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

Package: Single Pair

Size	Product Code
Small	592RPRS
Medium	592RPRM
Large	592RPRL
X-Large	592RPRX

Donning Aids



Sim-Slide®



Video: How to Use the Sim-Slide
To watch the video go to: bit.ly/SIGssyt

Sim-Slide

Size	Product Code
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

Sizing Chart

For use with Open toe styles

Women's Shoe Size	≤4–6	6.5–8	8.5–9.5	≥10
Men's Shoe Size	≤3–4.5	5–6.5	7–8	≥8.5
Size	S	M	L	XL



Magnide On/Off®



Video: How to Use the Magnide On/Off
To watch the video go to: bit.ly/SIGmofyt

Magnide On/Off

Size	Product Code
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

Sizing Chart

For use with Closed toe styles

Women's Shoe Size	≤8.5	9–12	≥12.5
Men's Shoe Size	≤7.5	8–11	≥11.5
Size	M	L	XL*

* XL suitable for extra-wide feet

Donning Aid



Cotton Underliner

10–15mmHg stocking. Can be used as an underlayer to protect fragile skin, secure wound dressings, to add 10–15mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings. For sizing, please refer to the Natural Rubber Calf sizing on page 32.

Package: 2 per pack

Size	Product Code
Small / Short	588USS
Small / Long	588USL
Medium / Short	588UMS
Medium / Long	588UML
Large-XL / Short	588ULS
Large-XL / Long	588ULL

Washing Solution



Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

Package Options

Contains	Product Code
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975
One 2 oz. Bottle (60ml)	586W985

Body Adhesive



IT STAYS!®

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin. IT STAYS! is pliable and moves with the body.

Package Options

Contains	Product Code
Twelve 2 oz. Bottles (60ml)	581I300
One 2 oz. Bottle (60ml)	581I350



ACCESSORIES

Donning Aids



Foot Slip

Comfortable applicator for use with open toe stockings.

Package: 6 per pack

Size	Product Code
One Size	580F200



Sigvaris DOFF N' DONNER®

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Non-latex.

Size	Product Code
One Size	583D100



Sigvaris DOFF N' DONNER CONE

This loading device makes it easy to roll stockings onto the Sigvaris DOFF N' DONNER. The bottom of the base firmly attaches to most clean flat surfaces via suction.



Video: How to use the Sigvaris DOFF N' DONNER
To watch the video go to: bit.ly/sigdnd

Size	Product Code
One Size	583D200

Caution: this product contains natural rubber latex which may cause allergic reactions



MELANY®

The donning and doffing device for compression stockings in case of limited mobility.

Sizing Chart				
Size	X-SMALL	SMALL	MEDIUM	LARGE
Calf Circumference	10"-12" 25-30cm	12"-15" 30-38cm	15"-18" 38-46cm	18"-22" 46-56cm
Product Code	579MXS	579MSM	579MMD	579MLG

Measure at the widest part of the calf



Video: How to Use the MELANY
To watch the video go to: goto.sigvaris.com/bbiMdz

⚠ WARNING: This product can expose you to chemicals, including Carbon black, which are known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

GLOSSARY OF TERMS

Term	Definition	Sigvaris Usage
Chronic Venous Disease (CVD)	A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Insufficiency (CVI), skin issues, and venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures Medical and Well Being products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
Double-Covered Yarn	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
Deep Vein Thrombosis (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
Graduated Compression v. Anti-Embolism	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for bedridden patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
Hypoallergenic	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers non-latex products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO Certified	The International Organization for Standardization, a worldwide federation of national standards bodies from over 130 countries. The United States member body is identified as ANSI (American National Standards Institute).	SIGVARIS GROUP follows ISO Standard 13485:2016 for our graduated compression products manufactured in North and South America. For our MCE products manufactured in Holland, Michigan, we follow ISO 9001.
Microfiber	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural-fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in 840 Soft Opaque, 830, 143 & 183 Microfiber Patterns, and 820 Microfiber.
mmHg (Millimeters of Mercury)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
Merino Wool	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in 422 Thermoregulating Wool, 421 Merino Outdoor Socks, and 152 & 192 All-Season Merino Wool.
Sea Island Cotton	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our 220 Sea Island and 151 & 191 Sea Island Cotton Socks.
Sensinnov	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in 780 Sheer, 840 Soft Opaque, and 820 Microfiber to help keep thigh stockings in place all day long.
Supima Cotton	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS GROUP uses this premium cotton yarn in the 230 Cotton, 220 Sea Island Cotton and 151 & 191 Sea Island Cotton Socks, 146 & 186 Casual Cotton, and 360 Cushioned Cotton products, making them comfortable and ideal for those with sensitive skin.

SYMPTOMS & INDICATIONS

Symptoms / Indications	Compression Level at the Ankle
<ul style="list-style-type: none"> • Tired, achy legs • Occupational or evening edema • Leg discomfort from long hours of standing or sitting • Prevent varicosities & swelling during pregnancy • Preventative for legs predisposed to risk • Reduce swelling during travel • Preventative for DVT during travel 	15–20mmHg
<ul style="list-style-type: none"> • Patient with diabetes who experiences edema (see Cautions below) 	18–25mmHg
<ul style="list-style-type: none"> • Heavy, fatigued, aching legs • Mild edema in lower extremities • Mild varicosities with minimal edema • Mild varicosities during pregnancy • Lymphedema upper extremities • Post-procedure of small veins • Preventative for or management of DVT 	20–30mmHg
<ul style="list-style-type: none"> • Moderate varicosities with mild to moderate edema (also during pregnancy) • Post-procedure of larger veins to maintain treatment success • CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active and / or recurrent ulcerations) • Management and treatment of DVT or Post-Thrombotic Syndrome • Superficial phlebitis • Lymphedema after decongestant therapy to maintain reduction • Orthostatic / postural hypotension • Post-fracture, post-traumatic edema <p><i>Note: A prescription is recommended for stockings with the indications of 30mmHg or higher</i></p>	30–40mmHg
<ul style="list-style-type: none"> • Lymphedema after decongestant therapy to maintain reduction • Post-Thrombotic Syndrome • Varicosities and / or edema 	40mmHg or higher
<ul style="list-style-type: none"> • Severe arterial insufficiency, intermittent claudication, ischemia • Uncontrolled congestive heart failure 	Contraindications
<ul style="list-style-type: none"> • Skin sensitivities or allergies • Advanced neuropathy with sensory loss in leg or foot • Diabetes with arterial compromise • Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician • Eczema, hypodermatitis, dermatosclerosis, dermatitis 	Cautions
<p><i>No liability accepted for non-observance of contraindications and cautions.</i></p> <p>Rabe E, Földi E, et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC) 2021. Recommendations 29 and 30. pg S46.</p> <p>Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151.</p> <p>Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p> <p>Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery.</p>	

DEALER PORTAL

The SIGVARIS GROUP Dealer Portal is your online location for securely ordering SIGVARIS GROUP products and managing your account. Place and track your orders, review order history, chat with Customer Care, and much more at sigvaris-online.com.

Dealer Portal benefits

- 24/7 access-order any time!
- Sizing Wizard
- Three ways to securely place an order online: Catalog, Quick Entry or Excel Upload
- Access to order history and tracking
- Shipping options including drop ship and easy access to an address table for quick shipment of reorders
- Check your account balance and print invoices
- Credit card payments
- New product launch information
- Current promotions
- Inelastic compression order forms

Email sigvariscare@sigvaris.com and you will receive your active account within one business day!

Can't find the product you're searching for?

Type the product SKU in the search function and select for ordering.



EDUCATION

SIGVARIS GROUP has several new and exciting educational opportunities for 2024 to help you get a leg up on leg health knowledge. These opportunities include online learning, webinars, and custom training sessions to educate fitters, clinicians, and therapists. Additionally, we offer a variety of courses with opportunities to earn continuing education credits.

Learn more at sigvariseducation.com



ORDERING INSTRUCTIONS

How to order Medical products					
Series Number 23 = Cotton (See individual Catalog pages for series number)	Compression Level 1 = 15–20mmHg 2 = 20–30mmHg* 3 = 30–40mmHg 4 = 40–50mmHg** 5 = 50–60mmHg**	Product Style C = Knee-high or Calf (See individual Catalog pages for available styles)	Size ML = Medium Long (See individual Catalog pages for sizing)	Gender M = Men W = Women O = Open toe	Color 99 = Black (See individual Catalog pages for available colors)
23	2	C	ML	M	99
<p>* 18–25mmHg available in Diabetic Compression Socks only</p> <p>** 40–50mmHg available in Secure stocking and Natural Rubber; 50–60mmHg available in Natural Rubber only</p>					
<p>To order Plus Calf with grip-top or Plus Size Pantyhose: add “/P” after the color in the ordering code (860P is available for women only)</p>					
<p>To order Calf Style with grip-top: add “/S” after the color in the ordering code (see product Catalog pages for availability)</p>					
<p>To order Thigh with Waist Attachment: add “/L” to indicate the Left leg or “/R” to indicate the Right leg after the color in the ordering code (860W and 500W only)</p>					
<p>To order Motion High Tech, Motion Thermoregulating Wool, Armsleeves, Gloves, and Gauntlets: the letter for Gender is not used in the ordering code. <i>Example: 412CSM59 = Motion High Tech, 20–30mmHg, Closed toe Calf, Size Small Medium, Color Red</i></p>					

How to order Well Being products			
Series Number 120 = Sheer Fashion (See individual Catalog pages for series number)	Product Style N = Thigh-high (See individual Catalog pages for available styles)	Size A = Small (See individual Catalog pages for sizing)	Color 33 = Honey (See individual Catalog pages for available colors)
120	N	A	33
<p>To order Athletic Recovery Socks, Performance Sleeves, Merino Outdoor Socks, and Eversoft Diabetic Socks: Size is only one letter and the letter for Gender is not used in the ordering code. <i>Example: 421CL31 = Merino Outdoor Socks, 15–20mmHg, Closed toe Calf, Size Large, Color Olive</i></p>			

To order:

Call: 1.800.322.7744, ext. 1
Fax: 1.800.481.5488

Online Dealer Portal:
sigvaris-online.com

MEASURING INSTRUCTIONS

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as soon after waking as possible when swelling is minimal.

Calf style Open/Closed Toe	
1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is Small, Medium, Large, or Extra Large. <i>Note: This should be your first measurement.</i>
2	Measure the calf circumference at the fullest part of the calf.
3	Determine the calf length by measuring from the floor up to the fibular head. This measurement is taken on the outside of the leg.

Pantyhose And Thigh-High style Open/Closed Toe	
1	Measure the ankle circumference directly above the ankle bone. Use the appropriate series sizing chart to determine if ankle is Small, Medium, Large, or Extra Large. <i>Note: This should be your first measurement.</i>
2	Measure the calf circumference at the fullest part of the calf.
3	Skip Step 3 (but refer to General Guidelines #6 below).
4	Measure the thigh at its widest circumference.
5	Determine the leg length by measuring from the floor up to the gluteal fold.
6	Determine the hip circumference by measuring the widest part of the hips.

General guidelines	
1	Measurements should always be taken as soon after waking as possible, when swelling is minimal.
2	Measure against bare skin, not over clothing.
3	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to the right for examples).
4	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
5	Length measurements should be taken from the ground up.
6	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.



1 Ankle circumference directly above ankle bone



2 Calf circumference at fullest part of the calf



3 Calf length from the floor to the fibular head



4 Widest circumference of the thigh



5 Leg length from the floor up to the gluteal fold



6 Circumference at the hip

SIGVARIS
GROUP

Ordering:

Online

sigvaris-online.com

Call

1-800-322-7744, ext. 1

Fax

1-800-481-5488

Customer Care

Monday–Friday

8am–7pm (EST)

Visit us on the web:

sigvaris.com



SIGVARIS, Inc.

1119 Highway 74 South
Peachtree City, GA 30269

Tel 800-322-7744



SIGVARIS GROUP USA